


Pe Rabbit Anti-Human Cd28 Monoclonal Antibody For Flow Cytometry - P10747

Item Number: CM0007950



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD28 (Tp44), a key T-cell costimulatory receptor. Suitable for flow cytometry applications in human samples to study T-cell activation, proliferation, and cytokine production.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD28 mAb
Remarks/Alias	Tp44
Species	Human
Gene ID (Human)	940
Gene ID	940
Immunogen	Recombinant protein of human CD28 Full length (P10747).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P10747
Protein Weight	25kDa
Shipping	Ice bag