


Pe Rabbit Anti-Human Cd282/Tlr2 Monoclonal Antibody For Flow Cytometry - O60603

Item Number: CM0008074



Introduction

PE-conjugated rabbit monoclonal antibody targeting human TLR2/CD282 for flow cytometry (FC). Validated for human cross reactivity, ideal for innate immunity and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD282/TLR2 mAb
Remarks/Alias	TIL4; CD282
Species	Human
GeneID (Human)	7097
GeneID	7097
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 290-390 of human TLR2 (NP_003255.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O60603
Protein Weight	90kDa
Shipping	Ice bag