

Pe Rabbit Anti-Human Cd300 A/C Mab For Fc - Q9Ugn4 / Q08708

Item Number: CM0012054



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD300a/c (IRC1/IRC2, CLM-8) for flow cytometry. Validated for human cross-reactivity. Suitable for immunophenotyping and functional studies of immune cells.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD300 a/c mAb
Remarks/Alias	IRC1; IRC2; CLM-8; IRp60; IGSF12; CMRF35H; CMRF-35H; CMRF35-H; CMRF35H9; CMRF35-H9; IRC1/IRC2; CMRF-35-H9
Species	Human
GeneID (Human)	11314/10871
GeneID	11314/10871
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-178 of human CD300a (NP_009192.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9UGN4/Q08708
Protein Weight	33kDa
Shipping	Ice bag