

## Pe Rabbit Anti-Human Cd301/Clec10A Monoclonal Antibody For Fc - Q8Iun9

Item Number: CM0011591



### Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD301/CLEC10A, a C-type lectin receptor, for flow cytometry. Validated in human specimens, enabling detection of myeloid cells and immune profiling.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | PE Rabbit anti-Human CD301/CLEC10A mAb                                                                                    |
| Remarks/Alias    | HML; MGL; HML2; CD301; CLECSF13; CLECSF14                                                                                 |
| Species          | Human                                                                                                                     |
| Gene ID (Human)  | 10462                                                                                                                     |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 61-316 of human CD301/CLEC10A (NP_878910.1) |
| Source           | Rabbit                                                                                                                    |
| Category         | Monoclonal Antibodies                                                                                                     |
| Application      | FC                                                                                                                        |
| Cross Reactivity | Human                                                                                                                     |
| SWISS            | Q8IUN9                                                                                                                    |
| Protein Weight   | 28kDa/32kDa/35kDa                                                                                                         |
| Shipping         | Ice bag                                                                                                                   |