

Pe Rabbit Anti-Human Cd303/Bdca-2 Monoclonal Antibody For Fc - Q8Wtt0

Item Number: CM0011006



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD303/BDCA-2 (CLEC4C), validated for flow cytometry. This antibody recognizes a C-type lectin receptor on plasmacytoid dendritic cells involved in antigen capture and interferon regulation.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD303/BDCA-2 mAb
Remarks/Alias	DLEC; HECL; BDCA2; CD303; BDCA-2; CLECSF7; CLECSF11; PRO34150
Species	Human
GeneID (human)	170482
GeneID	170482
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 45-213 of human CD303/BDCA-2 (NP_569708.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8WTT0
Protein Weight	21kDa/25kDa
Shipping	Ice bag