


Pe Rabbit Anti-Human Cd32 Mab - P12318

Item Number: CM0011272



Introduction

PE-conjugated rabbit anti-human CD32 monoclonal antibody for flow cytometry (FC). Recognizes the low affinity IgG Fc receptor FCGR2A (CD32A, FcγRIIa) on myeloid cells. Suitable for human samples.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD32 mAb
Remarks/Alias	CD32; FCG2; FcGR; CD32A; CDw32; FCGR2; IGFR2; FCGR2A1; FcγgammaRIIa
Species	Human
GeneID (Human)	2212
GeneID	2212
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 36-218 of human CD32 (NP_001129691.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P12318
Protein Weight	35kDa
Shipping	Ice bag