

Pe Rabbit Anti-Human Cd33 Monoclonal Antibody For Fc - P20138

Item Number: CM0009428



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD33 (Siglec-3) for flow cytometry. Reacts with human monocytic and myeloid lineage cells. Ideal for immunology and immuno-oncology research. – P20138

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD33 mAb
Remarks/Alias	p67; SIGLEC3; SIGLEC-3
Species	Human
GeneID (Human)	945
GeneID	945
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-259 of human CD33 (NP_001763.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P20138
Protein Weight	40kDa
Shipping	Ice bag