

Pe Rabbit Anti-Human Cd335/Nkp46 Monoclonal Antibody For Fc - O76036

Item Number: CM0010251



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD335/NKp46 (Natural cytotoxicity triggering receptor 1). Validated for flow cytometry, ideal for NK cell phenotyping, activation studies, and tumor immunology research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD335/NKp46 mAb
Remarks/Alias	LY94; CD335; NKP46; NK-p46
Species	Human
Gene ID (Human)	9437
Gene ID	9437
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O76036
Protein Weight	21kDa/22kDa/23kDa/32kDa/34kDa
Shipping	Ice bag