

Pe Rabbit Anti-Human Cd336/Nk-P44 Mab For Fc - O95944

Item Number: CM0009208



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD336/NK-p44, validated for flow cytometry (FC). Recognizes the extracellular domain of this NK cell-activating receptor; useful for tumor immunology research. Expected molecular weight: 31 kDa.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD336/NK-p44 mAb
Remarks/Alias	LY95; CD336; NKP44; NK-p44; dJ149M18.1
Species	Human
GeneID (human)	9436
GeneID	9436
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-190 of human CD336/NK-p44 (NP_004819.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O95944
Protein Weight	31kDa
Shipping	Ice bag