

Pe Rabbit Anti-Human Cd34 Monoclonal Antibody For Flow Cytometry - P28906

Item Number: CM0008111

Introduction

PE Rabbit anti-Human CD34 Monoclonal Antibody for flow cytometry. Detects human CD34, a hematopoietic progenitor cell marker, on stem cells and endothelial cells. Suitable for FC applications. SWISS: P28906.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD34 mAb
Remarks/Alias	CD34; CD34 molecule; GIG3; MORT1
Species	Human
GeneID	947 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-290 of human CD34 (NP_001020280.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28906
Protein Weight	35kDa/41kDa
Shipping	Ice bag