


Pe Rabbit Anti-Human Cd352/Slamf6 Monoclonal Antibody For Fc - Q96Du3

Item Number: CM0011073



Introduction

PE-conjugated rabbit monoclonal antibody specific for human CD352/SLAMF6. Designed for flow cytometry applications in immune cell research. No cross-reactivity with other species reported. 24-37 kDa.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD352/SLAMF6 Monoclonal Antibody
Remarks/Alias	KALI, NTBA, CD352, KALib, Ly108, NTB-A, SF2000
Species	Human
Gene ID (Human)	114836
Gene ID	114836
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-215 of human CD352/SLAMF6 (NP_001171643.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q96DU3
Protein Weight	24kDa/37kDa
Shipping	Ice bag