

Pe Rabbit Anti-Human Cd361 Monoclonal Antibody For Fc - P34910

Item Number: CM0009415

Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD361 (EVI2B, EVDB) for flow cytometry. Cross-reacts with human samples; suitable for granulocyte differentiation and hematopoietic progenitor cell studies.

[Learn More](#)

Parameter	Specification
Product Name	PE Rabbit anti-Human CD361 mAb
Remarks/Alias	EVDB; CD361; D17S376
Species	Human
Gene ID (Human)	2124
Immunogen	Recombinant fusion protein containing amino acids 22-199 of human CD361 (NP_006486.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P34910
Protein Weight	49kDa
Shipping	Ice bag