


Pe Rabbit Anti-Human Cd38 Monoclonal Antibody For Fc - P28907

Item Number: CM0008369



Introduction

PE-conjugated rabbit monoclonal antibody against human CD38, designed for flow cytometry. Specifically recognizes human CD38, a multifunctional enzyme and receptor involved in calcium signaling and immune regulation. Ideal for immunological research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD38 mAb
Remarks/Alias	ADPRC1; cADPR1; ADPRC 1
Species	Human
Gene ID (Human)	952
Gene ID	952
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-300 of human CD38 (P28907).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28907
Protein Weight	34kDa
Shipping	Ice bag