

Pe Rabbit Anti-Human Cd39/Entpd1 Monoclonal Antibody For Flow Cytometry - P49961

Item Number: CM0008911

Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD39/ENTPD1 (CD39, NTPDase1). Validated for flow cytometry (FC). Ideal for studying ectonucleotidase function, immune regulation, and platelet aggregation.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD39/ENTPD1 mAb
Remarks/Alias	CD39; SPG64; ATPDase; NTPDase-1
Species	Human
GeneID (Human)	953
GeneID	953
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 255-442 of human CD39/ENTPD1 (NP_001767.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P49961
Protein Weight	34kDa/42kDa/46kDa/57kDa/58kDa/59kDa
Shipping	Ice bag