

Pe Rabbit Anti-Human Cd4 Monoclonal Antibody For Fc - P01730

Item Number: CM0009094



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD4 (T4/Leu-3), validated for flow cytometry. Essential for immune cell profiling, T-helper cell research, and HIV studies. Swiss-Prot: P01730.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD4 mAb
Remarks/Alias	T4; IMD79; Leu-3; OKT4D; CD4mut
Species	Human
GeneID (Human)	920
GeneID	920
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-207 of human CD4 (NP_000607.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01730
Protein Weight	51kDa
Shipping	Ice bag