

Pe Rabbit Anti-Human Cd42B Monoclonal Antibody For Flow Cytometry - P07359

Item Number: CM0010136



Introduction

PE-conjugated rabbit monoclonal antibody against human CD42b (GP1BA) for flow cytometry. Suitable for platelet biology and hemostasis research. Recognizes human CD42b; protein weight: 71 kDa.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD42b mAb
Remarks/Alias	GP1BA; BDPLT1; BDPLT3; BSS; CD42B; CD42b-alpha; DBPLT3; GP1B; GPIbA; GPIbalpha; VWDP; platelet glycoprotein Ib alpha chain
Species	Human
GeneID (Human)	2811
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P07359
Protein Weight	71kDa
Shipping	Ice bag