

Pe Rabbit Anti-Human Cd42D/Gpv Mab For Fc - P40197

Item Number: CM0008463

Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD42d/GPV (P40197), validated for flow cytometry (FC). Part of the GPIIb-V-IX complex mediating platelet adhesion and hemostasis.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD42d/GPV mAb
Remarks/Alias	GPV; CD42d
Species	Human
GeneID (Human)	2814
GeneID	2814
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 17-477 of human CD42d (NP_004479.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P40197
Protein Weight	61 kDa
Shipping	Ice bag