

Pe Rabbit Anti-Human Cd43 Monoclonal Antibody For Fc - P16150

Item Number: CM0005033

Introduction

PE-conjugated rabbit anti-human CD43 monoclonal antibody validated for flow cytometry (FC). Recognizes human leukosialin, a major leukocyte surface sialoprotein involved in T-cell activation and migration. Ideal for immunology and lymphocyte profiling studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD43 mAb
Remarks/Alias	LSN; CD43; GALGP; GPL115; LEU-22
Species	Human
GeneID (human)	6693
GeneID	6693
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human CD43 (NP_003114.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P16150
Protein Weight	40 kDa
Shipping	Ice bag