


Pe Rabbit Anti-Human Cd44 Monoclonal Antibody For Fc - P16070

Item Number: CM0010630



Introduction

PE-conjugated rabbit monoclonal antibody against human CD44, suitable for flow cytometry (FC). Recognizes human CD44 antigen; ideal for studying cell adhesion, migration, and lymphocyte homing in immunology and cancer research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD44 mAb
Remarks/Alias	CD44; CDW44; CSPG8; ECMR-III; HCELL; HUTCH-I; IN; LHR; MC56; MDU2; MDU3; MIC4; Pgp1; CD44 antigen
Species	Human
Gene ID (Human)	960
Gene ID	960
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-178 of human CD44 (NP_000601.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P16070
Protein Weight	82kDa
Shipping	Ice bag