

Pe Rabbit Anti-Human Cd44V6 Monoclonal Antibody For Fc - P16070

Item Number: CM0012028



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD44V6 for flow cytometry. Detects CD44 antigen, a hyaluronate receptor involved in cell adhesion, lymphocyte homing, and tumor progression. Ideal for immunology and oncology research applications.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD44V6 mAb
Remarks/Alias	IN; LHR; MC56; MDU2; MDU3; MIC4; Pgp1; CDW44; CSPG8; H-CAM; HCELL; ECM-III; HUTCH-1; HUTCH-I; ECMR-III; Hermes-1
Species	Human
GeneID (Human)	960
GeneID	960
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 223-267 of human CD44V6 (P16070)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P16070
Protein Weight	82kDa
Shipping	Ice bag