

Pe Rabbit Anti-Human Cd48 Monoclonal Antibody For Flow Cytometry - P09326

Item Number: CM0010296



Introduction

PE-conjugated rabbit monoclonal antibody against human CD48, a GPI-anchored immune regulator, designed for flow cytometry. Ideal for studying lymphocyte activation, NK cell cytotoxicity, and immunological synapse formation.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD48 mAb
Remarks/Alias	BCM1; BLAST; hCD48; mCD48; BLAST1; SLAMF2; MEM-102
Species	Human
Gene ID	962
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-220 of human CD48 (NP_001769.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P09326
Protein Weight	19kDa/28kDa
Shipping	Ice bag