


Pe Rabbit Anti-Human Cd5 Mab For Fc - P06127

Item Number: CM0004995



Introduction

PE Rabbit anti-Human CD5 monoclonal antibody for flow cytometry. Recognizes human CD5 (T1/Leu-1), a lymphoid-specific receptor on T cells and B1a cells. Suitable for immunophenotyping and immune regulation studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD5 mAb
Remarks/Alias	CD5; LEU1; T1; CD5 molecule
Species	Human
GeneID (Human)	921
GeneID	921
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-371 of human CD5 (NP_055022.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06127
Protein Weight	55kDa
Shipping	Ice bag