


Pe Rabbit Anti-Human Cd52 Monoclonal Antibody For Flow Cytometry - P31358

Item Number: CM0012305



Introduction

PE-conjugated rabbit monoclonal antibody against human CD52 (CAMPATH-1, HE5) for flow cytometry. Specific to human CD52, a GPI-anchored glycoprotein (~7 kDa) on lymphocytes and epididymal cells. Suitable for immunology and reproductive biology research.

[Learn More](#)

Specification	Details
Product Name	PE Rabbit anti-Human CD52 mAb
Remarks/Alias	HE5; CDW52; EDDM5
Species	Human
GeneID (Human)	1043
GeneID	1043
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-61 of human CD52 (NP_001794.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P31358
Protein Weight	7kDa
Shipping	Ice bag