


Pe Rabbit Anti-Human Cd55 Monoclonal Antibody For Fc - P08174

Item Number: CM0005328



Introduction

PE-conjugated rabbit monoclonal antibody against human CD55 (DAF, CR) for flow cytometry. Recognizes human CD55, a complement regulator and virus receptor. Suitable for human samples. UniProt P08174.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD55 mAb
Remarks/Alias	CR; TC; DAF; CROM; CHAPLE
Species	Human
Gene ID (Human)	1604
Gene ID	1604
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08174
Protein Weight	39kDa/41kDa/48kDa/56kDa/59kDa
Shipping	Ice bag