

Pe Rabbit Anti-Human Cd6 Monoclonal Antibody For Flow Cytometry - P30203

Item Number: CM0010448



Introduction

PE-conjugated rabbit monoclonal antibody against human CD6, for flow cytometry. Immunogen: recombinant human CD6 (aa 18-364). Reacts with human T cells, NK cells, B-1 cells. Research use only.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD6 mAb
Remarks/Alias	CD6; TP120; CD6 molecule
Species	Human
GeneID (Human)	923
GeneID	923
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-364 of human CD6 (NP_006716.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P30203
Protein Weight	58kDa/60kDa/63kDa/64kDa/68kDa/71kDa
Shipping	Ice bag