

Pe Rabbit Anti-Human Cd61 Monoclonal Antibody For Fc - P05106

Item Number: CM0008620



Introduction

PE Rabbit anti-Human CD61 Monoclonal Antibody for flow cytometry. Targets integrin beta-3 (GPIIIa), involved in platelet aggregation and hemostasis. Ideal for thrombosis research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD61 mAb
Remarks/Alias	GT; GT2; CD61; GP3A; BDPLT2; GPIIIa; BDPLT16; BDPLT24
Species	Human
GeneID (Human)	3690
GeneID	3690
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-718 of human CD61 (NP_000203.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P05106
Protein Weight	86kDa/87kDa
Shipping	Ice bag