


Pe Rabbit Anti-Human Cd62E/E-Selectin Mab For Fc - P16581

Item Number: CM0012771



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD62E/E-Selectin for flow cytometry. Validated for immunoadhesion and inflammation research, reacts specifically with human samples.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD62E/E-Selectin mAb
Remarks/Alias	ELAM; ESEL; CD62E; ELAM1; LECAM2; selectin-e
Species	Human
GeneID (Human)	6401
GeneID	6401
Immunogen	Recombinant protein, Recombinant fusion protein containing a sequence corresponding to amino acids 22-230 of human CD62E/E-Selectin (NP_000441.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P16581
Protein Weight	66kDa
Shipping	Ice bag