

Pe Rabbit Anti-Human Cd62L/L-Selectin Monoclonal Antibody For Fc - P14151

Item Number: CM0008864



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD62L/L-Selectin (SELL). Validated for flow cytometry (FC). Detects a 42kDa/43kDa protein involved in leukocyte adhesion and homing.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD62L/L-Selectin mAb
Remarks/Alias	SELL; CD62L; LAM1; LECAM1; LEU8; LNHHR; LSEL; LYAM1; PLNHR; TQ1; L-selectin
Species	human
Gene ID (Human)	6402
Gene ID	6402
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 39-332 of human CD62L/L-Selectin (NP_000646.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P14151
Protein Weight	42kDa/43kDa
Shipping	Ice bag