

Pe Rabbit Anti-Human CD62P/P-Selectin Monoclonal Antibody For Fc - P16109

Item Number: CM0010706



Introduction

This PE-conjugated rabbit monoclonal antibody targets human CD62P (P-selectin), a leukocyte adhesion receptor involved in inflammation. Validated for flow cytometry (FC) with human cross-reactivity. Suitable for immunological and inflammation research studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD62P/P-selectin mAb
Remarks/Alias	CD62; GRMP; PSEL; CD62P; GMP140; LECAM3; PADGEM
Species	Human
GeneID	6403
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 42-771 of human CD62P/P-selectin (NP_002996.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P16109
Protein Weight	91kDa
Shipping	Ice bag