


Pe Rabbit Anti-Human Cd63 Monoclonal Antibody For Fc - P08962

Item Number: CM0004456



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD63 (LAMP-3, ME491), validated for flow cytometry. Suitable for research in cell adhesion, integrin signaling, melanoma, and tetraspanin biology.

[Learn More](#)

Parameter	Specification
Product Name	PE Rabbit anti-Human CD63 mAb
Remarks/Alias	MLA1; ME491; LAMP-3; OMA81H; TSPAN30
Species	Human
Gene ID (Human)	967
Gene ID	967
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 1-238 of human CD63, NP_001771.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Human
SWISS	P08962
Protein Weight	26kDa
Shipping	Ice bag