

Pe Rabbit Anti-Human Cd7 Monoclonal Antibody For Flow Cytometry - P09564

Item Number: CM0010572

Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD7 (GP40, TP41). Validated for flow cytometry, detects T-cell and NK cell marker. Ideal for immunophenotyping and leukemia research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD7 mAb
Remarks/Alias	GP40, TP41, Tp40, LEU-9
Species	Human
GeneID (Human)	924
GeneID	924
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-180 of human CD7 (NP_006128.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P09564
Protein Weight	25kDa
Shipping	Ice bag