

Pe Rabbit Anti-Human Cd72 Monoclonal Antibody For Flow Cytometry - P21854

Item Number: CM0009547

Introduction

High-quality PE-conjugated rabbit anti-human CD72 monoclonal antibody validated for flow cytometry (FC). Ideal for immunophenotyping B-cell subsets and studying BCR signaling. Specific to human CD72.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD72 mAb
Remarks/Alias	CD72; CD72b; LYB2; CD72 molecule
Species	Human
GeneID	971 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 117-359 of human CD72 (NP_001773.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P21854
Shipping	Ice bag