

Pe Rabbit Anti-Human Cd79B/Igβ Monoclonal Antibody For Fc - P40259

Item Number: CM0011291



Introduction

Rabbit monoclonal antibody specific for human CD79b/Igβ, conjugated to PE for flow cytometry. Detects B-cell antigen receptor complex-associated protein beta chain in human B-cells. Swiss-Prot P40259.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD79b/Igβ mAb
Remarks/Alias	CD79B; AGM6; B29; IGB; CD79b molecule
Species	Human
GeneID (Human)	974
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 29-159 of human CD79b/Igβ (NP_000617.1))
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P40259
Protein Weight	14kDa/26kDa
Shipping	Ice bag