


Pe Rabbit Anti-Human Cd8 Polymab® Monoclonal Antibody For Fc - P01732 / P10966

Item Number: CM0010076



Introduction

PE-conjugated rabbit anti-human CD8 PolymAb® monoclonal antibody for flow cytometry (FC). Detects CD8 alpha and beta chains (UniProt: P01732, P10966) with human cross-reactivity. Molecular weight: 26 kDa/24 kDa.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD8 PolymAb®
Remarks/Alias	CD8; p32; Leu2; IMD116; CD8alpha; LY3; P37; LEU2; LYT3; Ly-3; CD8B1; CD8beta.
Species	Human
GenelD	925/926
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P01732/P10966
Protein Weight	26 kDa/24 kDa
Shipping	Ice bag