

Pe Rabbit Anti-Human Cd80/B7-1 Monoclonal Antibody For Fc - P33681

Item Number: CM0005506



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD80 (B7-1), validated for flow cytometry. CD80 is a key costimulatory molecule on antigen-presenting cells critical for T-cell activation and immune checkpoint research.

[Learn More](#)

Parameter	Specification
Product Name	PE Rabbit anti-Human CD80/B7-1 mAb
Remarks/Alias	B7, BB1, B7-1, B7.1, LAB7, CD28LG, CD28LG1
Species	Human
GeneID (human)	941
GeneID	941
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-242 of human CD80/B7-1 (NP_005182.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P33681
Protein Weight	33kDa
Shipping	Ice bag