


Pe Rabbit Anti-Human Cd81 Monoclonal Antibody For Flow Cytometry - P60033

Item Number: CM0004457



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD81, validated for flow cytometry. Detects CD81 on immune cells; 26kDa molecular weight. Ideal IVD raw material for research and diagnostic assay development.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD81 mAb
Remarks/Alias	S5.7; CVID6; TAPA1; TSPAN28
Species	Human
Gene ID (Human)	975
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P60033
Protein Weight	26kDa
Shipping	Ice bag