

Pe Rabbit Anti-Human Cd82 Monoclonal Antibody For Fc - P27701

Item Number: CM0009088



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD82 (KAI1/TSPAN27), a tetraspanin metastasis suppressor. Validated for flow cytometry. Suitable for immunophenotyping of lymphoid cells and cancer research applications.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD82 mAb
Remarks/Alias	R2; 4F9; C33; IA4; ST6; GR15; KAI1; SAR2; TSPAN27
Species	Human
Gene ID (Human)	3732
Gene ID	3732
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-267 of human CD82 (NP_002222.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P27701
Protein Weight	27kDa/30kDa
Shipping	Ice bag