

Pe Rabbit Anti-Human Cd85D/E/I/J Monoclonal Antibody For Fc - Q8N423 / Q8N6C8 / O75019 / Q8Nhl6

Item Number: CM0011645



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD85d/e/i/j (LILRB2, LILRA3, LILRB1, LILRB4). Validated for flow cytometry, cross-reacts with human samples. Ideal for immune tolerance and checkpoint research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD85d/e/i/j mAb
Remarks/Alias	ILT4; LIR2; CD85D; ILT-4; LIR-2; MIR10; MIR-10
Species	Human
GeneID (human)	10288/11026/11024/10859
GeneID	10288/11026/11024/10859
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-461 of human CD85d (NP_001074447.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8N423/Q8N6C8/O75019/Q8NHL6
Protein Weight	65kDa
Shipping	Ice bag