

Pe Rabbit Anti-Human Cd85I/J Monoclonal Antibody For Fc - O75019 / Q8Nhl6

Item Number: CM0009808

Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD85i/j (LILRA1, LIR6) for flow cytometry. Validated for human reactivity, recombinant protein immunogen, ~53kDa. Ideal for immunology research and immune profiling.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD85i/j mAb
Remarks/Alias	LIR6; CD85I; LIR-6
Species	Human
Gene ID (Human)	11024/10859
Gene ID	11024/10859
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O75019/Q8NHL6
Protein Weight	53kDa
Shipping	Ice bag