


Pe Rabbit Anti-Human Cd85j/Ilt2 (Lilrb1) Monoclonal Antibody For Flow Cytometry - Q8Nhl6

Item Number: CM0005003



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD85j/ILT2 (LILRB1), an inhibitory MHC class I receptor. Validated for flow cytometry, suitable for immune cell profiling and natural killer cell research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD85j/ILT2 mAb
Remarks/Alias	LILRB1; CD85j; ILT-2; ILT2; LIR-1; LIR1; MIR-7; MIR7; PIR-B; PIRB; leukocyte immunoglobulin like receptor B1
Species	Human
Gene ID	10859
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-458 of human CD85j/ILT2 (NP_006660.4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q8NHL6
Protein Weight	49kDa/70kDa/71kDa
Shipping	Ice bag