


Pe Rabbit Anti-Human Cd87/Plaur Monoclonal Antibody For Fc - Q03405

Item Number: CM0010830



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD87/PLAUR. Validated for flow cytometry. CD87/PLAUR is a GPI-anchored receptor that recruits uPA and mediates extracellular matrix degradation and cell adhesion.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD87/PLAUR mAb
Remarks/Alias	PLAUR; CD87; U-PAR; UPAR; URKR; plasminogen activator, urokinase receptor
Species	Human
Gene ID (Human)	5329
Gene ID	5329
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-303 of human CD87/PLAUR (NP_002650.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q03405
Protein Weight	37kDa31kDa/32kDa/36kDa
Shipping	Ice bag