

Pe Rabbit Anti-Human Cd8B Monoclonal Antibody For Flow Cytometry - P10966

Item Number: CM0007946

Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD8b (SWISS P10966). Validated for flow cytometry. Detects CD8 beta chain, a coreceptor for MHC class I essential for T-cell activation and thymic selection. Suitable for T-cell immunophenotyping and diagnostics research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD8b mAb
Remarks/Alias	LY3; P37; LEU2; LYT3; Ly-3; CD8B1; CD8beta
Species	Human
GeneID (Human)	926
GeneID	926
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 22-170 of human CD8b (NP_004922.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P10966
Protein Weight	24kDa
Shipping	Ice bag