


Pe Rabbit Anti-Human Cd9 Monoclonal Antibody For Flow Cytometry - P21926

Item Number: CM0009271



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD9 (tetraspanin-29), validated for flow cytometry. CD9 regulates cell adhesion, platelet activation, and exosome biology. Suitable for immunophenotyping hematopoietic and epithelial cells.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD9 mAb
Remarks/Alias	MIC3; MRP-1; BTCC-1; DRAP-27; TSPAN29; TSPAN-29
Species	Human
GeneID (Human)	928
GeneID	928
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-228 of human CD9 (NP_001760.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P21926
Protein Weight	25kDa
Shipping	Ice bag