

Pe Rabbit Anti-Human Cd90/Thy1 Mab - P04216

Item Number: CM0010091

Introduction

PE-conjugated rabbit anti-human CD90/Thy1 monoclonal antibody for flow cytometry. Targets human CD90 (Thy-1 membrane glycoprotein, UniProt P04216), a GPI-anchored cell surface protein implicated in cell-cell interactions and synaptogenesis.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD90/Thy1 mAb
Remarks/Alias	CD90; CDw90
Species	Human
GeneID (Human)	7070
GeneID	7070
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 20-130 of human CD90/Thy1 (NP_006279.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P04216
Protein Weight	18kDa
Shipping	Ice bag