

Pe Rabbit Anti-Human Cd94 Mab For Fc - Q13241

Item Number: CM0010642



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD94 (KLRD1) for flow cytometry. Enables detection of the NK cell receptor involved in immune surveillance, tumor immunity, and viral escape mechanisms in human samples.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD94 mAb
Remarks/Alias	CD94
Species	Human
GeneID (Human)	3824
GeneID	3824
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-179 of human CD94. (NP_002253.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q13241
Protein Weight	17kDa/20kDa
Shipping	Ice bag