

Pe Rabbit Anti-Human Cd96/Tactile Mab For Fc - P40200

Item Number: CM0009284



Introduction

PE-conjugated rabbit monoclonal antibody against human CD96/TACTILE protein. Validated for flow cytometry (FC). Detects CD96, a T-cell and NK cell adhesion molecule. Ideal for immune cell phenotyping and immunological research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD96/TACTILE mAb
Remarks/Alias	TACTILE
Species	Human
Gene ID (Human)	10225
Gene ID	10225
Immunogen	Recombinant protein corresponding to a sequence within amino acids 25-519 of human CD96/TACTILE (NP_937839.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P40200
Protein Weight	66kDa
Shipping	Ice bag