


Pe Rabbit Anti-Human Cdkn2A/P16Ink4A Monoclonal Antibody For Fc (Intra) - P42771 / Q8N726

Item Number: CM0009952



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CDKN2A/p16INK4a for intracellular flow cytometry. Detects a key tumor suppressor and cell cycle inhibitor, supporting cancer research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CDKN2A/p16INK4a mAb
Remarks/Alias	ARF; MLM; P14; P16; P19; CMM2; INK4; MTS1; TP16; CDK4I; CDKN2; INK4A; MTS-1; P14ARF; P19ARF; P16INK4; P16INK4A; P16-INK4A
Species	Human
GeneID (human)	1029
GeneID	1029
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-156 of human CDKN2A/p16INK4a (NP_000068.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P42771 / Q8N726
Protein Weight	16 kDa
Shipping	Ice bag