

Pe Rabbit Anti-Human Cpt1A Monoclonal Antibody For Intracellular Flow Cytometry - P50416

Item Number: CM0008721



Introduction

PE-conjugated rabbit monoclonal antibody targeting human CPT1A (CPT1-L), validated for intracellular flow cytometry. Recombinant immunogen (aa 520-720); detects human CPT1A. Suitable for fatty acid metabolism and mitochondrial studies.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CPT1A mAb
Remarks/Alias	CPT1; CPT1-L; L-CPT1
Species	Human
GeneID (Human)	1374
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 520-720 of human CPT1A (NP_001867.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	Flow Cytometry (Intracellular)
Cross Reactivity	Human
SWISS	P50416
Protein Weight	88kDa
Shipping	Ice bag