

Pe Rabbit Anti-Human Desmoglein 3 (Dsg3) Monoclonal Antibody For Flow Cytometry - P32926

Item Number: CM0005547



Introduction

PE-conjugated rabbit monoclonal antibody targeting human Desmoglein 3 (DSG3/PVA/CDHF6) for flow cytometry. Validated for human cross-reactivity, this antibody aids cell adhesion and desmosome research.

[Learn More](#)

Parameter	Specification
Product Name	PE Rabbit anti-Human Desmoglein 3/PVA mAb
Remarks/Alias	PVA; ABOLM; CDHF6
Species	Human
GeneID (Human)	1830
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 50-615 of human Desmoglein 3/PVA, NP_001935.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32926
Protein Weight	108kDa
Shipping	Ice bag