


Pe Rabbit Anti-Human Egfr Mab For Flow Cytometry - P00533

Item Number: CM0011309



Introduction

PE-conjugated rabbit monoclonal antibody targeting human EGFR (ERBB1). Validated for flow cytometry (FC) to analyze EGFR expression in cell surface and intracellular signaling studies. Ideal for cancer and receptor tyrosine kinase research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human EGFR mAb
Remarks/Alias	ERBB; ERRP; HER1; mENA; ERBB1; PIG61; NISBD2
Species	Human
Gene ID (Human)	1956
Gene ID	1956
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 335-525 of human EGFR(NP_005219.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P00533
Protein Weight	44kDa/69kDa/77kDa/134kDa
Shipping	Ice bag